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REF 122 500

D12250012020-031002 ML

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EN BORT EpiContur with 2 Pads

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a brace for relief of the muscle origins on the ulnar and/or radial epicondyle.

Indications

Epicondylitis, chronic or recurring epicondylopathy, tendopathies, excessive strain, so-called tennis/golf arm.

Contraindications

Peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the area of the body treated.

Application risks/Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- remove it during longer periods of rest, e.g. when sleeping
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use – this medical aid is intended for treating one patient
- whilst wearing the brace, please neither use any creams nor ointments on or around it as they can destroy the material

Fitting instructions

The pads can be placed freely using the Velcro fastener. For correct positioning, observe the markings printed on the support. The narrow part of the pad should always point in the direction of the elbow joint. Firstly, thread the Velcro tape through the plastic eyelet and then slip the clasp over the lower arm. The pads should exert pressure on the muscle belly which leads from the hand to the elbow, approx. 5 cm beneath the elbow joint. The muscle belly can be clearly felt when you make a fist. The BORT epicondylitis clasp should only exert pressure on this muscle if you make a fist or grasping movement. When the hand is relaxed, the clasp should be placed on the pads without noticeable pressure. Close the Velcro fastener so that the clasp fits closely but causes no congestion. Cut off any protruding Velcro tape.

Removal

To remove the epicondylitis brace, loosen the Velcro tape, increase the circumference and remove the brace in the direction of the hand.

Material composition

Polyamide (PA), polyurethane (PUR)
The sewn in textile label provides the precise material composition.

Cleaning information

Hand-washing Do not clean chemically Do not bleach
 Do not dry in a tumble dryer Do not iron
Do not use fabric conditioner. Dry in the fresh air.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the appointed authority for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 11.2019

Medical device | Single patient – multiple use

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